

# BPX5 - A Capillary Column which will perform 80 % of all applications

## WHAT IS BPX5?

BPX5 is a 5 % phenyl (equivalent) / 95 % methyl polysilphenylene / siloxane phase.

## WHY USE A 5 % PHENYL PHASE?

It is estimated that 80 % of all analyses can be carried out on one column with one set of dimensions and film thickness. The 5 % phenyl column is still a non-polar column, but the small amount of phenyl content gives the phase a little edge in selectivity compared with the 100 % methyl phases.

The non-polar phases tend to have lower bleed and are more robust compared with polar phases.

## WHAT CAN THE BPX5 COLUMN BE USED FOR?

The BPX5 is a multi-purpose column which can be used for a range of applications from non-polar components (such as those in gasoline) to polar components such as amines, pesticides, triglycerides, drugs and polycyclic aromatic hydrocarbons.

## EQUIVALENT COLUMNS


If you are using any 5 % phenyl column, BPX5 can replace it. The table below is an excellent replacement reference.

This column is a **direct replacement** for;

CP-SIL 8CB	MDN-5
CP-SIL 8CBms	PE-5
DB-5	PTE-5
DB5.625	Rtx-5
DB5-MS	Rtx-5MS
HP-5	SPB-5
HP-5MS	XTI-5
HP5-TA	ZB-5 (RH-5)
MDN-5S	007-5
HP-Ultra2	007-5MS


## OPERATING TEMPERATURES

0.25 µm – 1.0 µm	-40 °C to 360 °C/370 °C
>1.0 µm	-40 °C to 350 °C/360 °C




Phenols


# PRODUCT DATA




Pesticides - broad range pesticide mixture




Analysis of chlorine, phosphorous and nitrogen based pesticides



Pesticides - EPA method 608 and additional pesticides




USEPA method 608.1 - organochlorine pesticides




USEPA method 614 - organophosphorous pesticides


# PRODUCT DATA




Triglycerides - analysis of mono, di, and triglycerides




Triglyceride distribution in milk fat



Standard mixture of triglycerides




Amines - analysis of aniline impurities




Analysis of aromatic amines


# PRODUCT DATA



Fragrance - analysis of lavender essential oil




Analysis of rosemary essential Oil



Analysis of eucalyptus essential Oil

Phase: BPX5, 1.0  $\mu$ m  
 Column: 25 m x 0.53 mm ID  
 Initial Temp.: 120  $^{\circ}$ C  
 Rate: 10  $^{\circ}$ C/min  
 Final Temp.: 310  $^{\circ}$ C  
 Detector: FID  
 Injector: Split, 240  $^{\circ}$ C  
 Carrier Gas: H<sub>2</sub>, 2 psi  
 Part No.: 054131




1. Methamphetamine
2. Phendimetrazine
3. Phencyclidine
4. Mepivocaine
5. Methaqualone
6. Amitriptyline
7. Codeine
8. Diazepam
9. Fentanyl

Note: A high degree of inertness and the high maximum operating temperature (for semi-volatile drugs), makes the BPX5 column suitable for this class of compounds.

Pharmaceuticals - basic drug screen (10-20 ng/component)

Phase: BPX5, 0.25  $\mu$ m  
 Column: 25 m x 0.22 mm ID  
 Initial Temp.: 180  $^{\circ}$ C  
 Rate: 8  $^{\circ}$ C/min  
 Final Temp.: 350  $^{\circ}$ C, 10 min  
 Detector: FID  
 Sensitivity: 32 x 10<sup>-12</sup> AFS  
 Injection Mode: Split  
 Carrier Gas: H<sub>2</sub>, 10 psi  
 Part No.: 054113




1. Testosterone
2. Pregnenolone
3. Progesterone
4. Cholesterol

Note: Whether your particular application is forensic screening or pharmaceutical testing, the BPX5 phase is ideal.

Steroid analysis (underivatised)

Phase: BPX5, 0.5  $\mu$ m film  
 Column: 25 m x 0.32 mm ID  
 Initial Temp.: 180  $^{\circ}$ C, 0 min  
 Rate: 15  $^{\circ}$ C/min  
 Final Temp.: 350  $^{\circ}$ C, 5 min  
 Detector: FID  
 Sensitivity: 140 x 10<sup>-12</sup> AFS  
 Injection Mode: Split  
 Carrier Gas: H<sub>2</sub>, 6 psi  
 Part No.: 054125




1. Dimethyl Phthalate
2. Diethyl Phthalate
3. Diisopropyl Phthalate
4. Dipropyl Phthalate
5. Diisobutyl Phthalate
6. Dibutyl Glutarate
7. Dipropyl Phthalate
8. Dihexyl Phthalate
9. Diheptyl Phthalate
10. Branched Dinonyl Phthalate
11. Dioctyl Phthalate
12. Diundecyl Phthalate
13. Didodecyl Phthalate

Phthalates - analysis of dimethyl esters of phthalate esters

# PRODUCT DATA

Phase: BPX5, 0.25 µm film  
 Column: 25 m x 0.22 mm ID  
 Initial Temp: -20 °C, 1 min  
 Rate 1: 4 °C/min  
 Temp 2: 40 °C  
 Rate 2: 8 °C/min  
 Final Temp: 150 °C  
 Detector: FID  
 Sensitivity: 10 x 10<sup>-12</sup> AFS  
 Injection Mode: Split  
 Carrier Gas: H<sub>2</sub>, 10 psi  
 Part No.: 054113

- |                            |                             |
|----------------------------|-----------------------------|
| 1. iso-Butane              | 17. 2,3-Dimethylhexane      |
| 2. n-Butane                | 18. Ethylbenzene            |
| 3. iso-Pentane             | 19. m-Xylene                |
| 4. n-Pentane               | 20. p-Xylene                |
| 5. 2,2-Dimethylbutane      | 21. o-Xylene                |
| 6. 2,3-Dimethylbutane      | 22. n-Propylbenzene         |
| 7. 2-Methylpentane         | 23. 1-Methyl-3-ethylbenzene |
| 8. 3-Methylpentane         | 24. 1-Methyl-4-ethylbenzene |
| 9. n-Hexane                | 25. 1,3,5- Trimethylbenzene |
| 10. Benzene                | 26. 1-Methyl-2-ethylbenzene |
| 11. 2-Methylhexane         | 27. 1,2,4- Trimethylbenzene |
| 12. 2,3-Dimethylpentane    | 28. 1,2,3-Trimethylbenzene  |
| 13. 3-Methylhexane         | 29. Naphthalene             |
| 14. 2,2,4-Trimethylpentane | 30. 2-Methylnaphthalene     |
| 15. n-Heptane              | 31. 1-Methylnaphthalene     |
| 16. Toluene                |                             |



Petroleum - leaded gasoline

## ORDERING INFORMATION

ID mm	Film µm	10 metre	15 metre	30 metre	60 metre
0.1	0.10	054099	-	-	-
0.22	0.25	-	-	054142	-
0.25	0.25	-	054100	054101	054102
0.25	0.50	-	-	0541025	-
0.25	1.00	-	054121	054122	054123
0.32	0.25	-	054144	054145	054146
0.32	0.50	-	-	0541205	-
0.32	1.00	-	054152	054153	054154
0.53	0.50	-	-	0541345	-
0.53	1.00	-	054147	054148	054158
0.53	1.50	-	0541347	0541348	-
0.53	3.00	-	054159	054160	-

ID mm	Film µm	12 metre	25 metre	50 metre
0.15	0.25	054103	054104	054105
0.15	0.40	054107	054108	-
0.22	0.25	054112	054113	054114
0.22	1.00	-	054116	054117
0.32	0.25	054118	054119	054120
0.32	0.50	054124	054125	054126
0.32	1.00	054127	054128	054129
0.32	3.00	-	054136	-
0.53	0.25	054133	054134	-
0.53	1.00	054130	054131	054132
0.53	3.00	054138	054139	-

\*Mentioned trademarks are registered trademarks of their respective owners.  
 Fused Silica capillary columns are manufactured under license granted by Agilent Technologies (HP).

For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

### AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd  
 Toll Free: 1800 800 167  
 Tel: +61 (0) 3 9837 4200  
 Fax: +61 (0) 3 9874 5672  
 Email: [support@sge.com](mailto:support@sge.com)

### CHINA

SGE Shanghai Representative Office  
 Tel: +86 21 6407 9382  
 Fax: +86 21 6407 9386  
 Email: [china@sge.com](mailto:china@sge.com)

### MIDDLE EAST

SGE Gulf  
 Tel: +971 6 557 3341  
 Fax: +971 6 557 3541  
 Email: [gulfsupport@sge.com](mailto:gulfsupport@sge.com)

### EUROPE

SGE Europe Ltd  
 European Head Office  
 Toll Free: 00800 2790 8999  
 Toll Free Fax: 00800 2626 2609  
 Tel: +44 1908 568 844  
 Fax: +44 1908 566 790  
 Tel France: +33 1 69 29 80 90  
 Fax France: +33 1 69 29 09 25  
 Tel Germany: +49 (0) 6155 / 60746 0  
 Fax Germany: +49 (0) 6155 / 60746 50  
 Email: [europe@sge.com](mailto:europe@sge.com)

### INDIA

SGE Laboratory Accessories Pvt Ltd  
 Tel: +91 22 24715896  
 Fax: +91 22 24716592  
 Email: [sgeindia@vsnl.com](mailto:sgeindia@vsnl.com)

### UNITED STATES OF AMERICA

SGE Incorporated  
 Toll Free: (800) 945 6154  
 Tel: +1 512 837 7190  
 Fax: +1 512 836 9159  
 Email: [usa@sge.com](mailto:usa@sge.com)

### JAPAN

SGE Japan Inc  
 Tel: +81 45 222 2885  
 Fax: +81 45 222 2887  
 Email: [japan@sge.com](mailto:japan@sge.com)

 **SGE Analytical Science**  
[www.sge.com](http://www.sge.com)



## BP5

ID (mm)	Film Thickness (µm)	Length (m)	Temperature Limits (°C)	Part No.
0.25	0.25	60	-60 to 320/340	054184
0.25	1	60	-60 to 320/340	054215
0.32	0.25	12	-60 to 320/340	054179
0.32	0.25	15	-60 to 320/340	054176
0.32	0.25	25	-60 to 320/340	054180
0.32	0.5	25	-60 to 320/340	054186
0.32	1	25	-60 to 320/340	054192
0.32	0.25	30	-60 to 320/340	054177
0.32	0.5	30	-60 to 320/340	054216
0.32	1	30	-60 to 320/340	054189
0.32	0.5	50	-60 to 320/340	054187
0.32	1	50	-60 to 320/340	054193
0.32	0.25	60	-60 to 320/340	054178
0.32	1	60	-60 to 320/340	054188
0.53	1	12	-60 to 320/340	054197
0.53	1	15	-60 to 320/340	054194
0.53	1.5	15	-60 to 320/340	054199
0.53	1	25	-60 to 320/340	054198
0.53	0.5	30	-60 to 320/340	0541935
0.53	1	30	-60 to 320/340	054195
0.53	5	30	-60 to 280/300	054196
0.53	1.5	60	-60 to 280/300	054204

### GC Columns and Applications

#### Expert Tip :

If the injection port temperature is not specified in the method, 250 °C is usually the recommended temperature.




## GC Capillary Columns | 5% Phenyl Polysilphenylene-siloxane

### BPX5

- High temperature.
- General purpose GC column – suitable for over 80% of all routine analyses performed by gas chromatography.
- Very low bleed – ideal for trace analysis.
- Non-polar.
- Extremely inert.
- Ideal for GC-MS.
- 360 – 370 °C upper temperature limit – dependent on film thickness.

Recommended column for General Purpose use.

Especially Suitable for these Industries:	 All Industries
Areas:	Ultra trace analyses, pesticides/herbicides, hydrocarbons, solvents, phenols, amines, GC/MS and other specific detector applications. Applications ENV62, ARO09, ENV20, ENV03, ENV48, ENV59, ENV84, FOO21, FLA14, FLA16, FLA15, FLA12 FLA13, ENV54, PET22, SOL33 PHA06, PHA08, PHA15.
Suitable Replacement for:	DB-5, DB-5ms, DB-5.625, XTI-5, Rtx-5ms, Ultra-2, HP-5, HP-5MS, HP5-TA, SPB-5, MDN-5S, CP-Sil8CB, Rxt-Sil 5MS, AT-5ms, VB-5, ZB-5, VF-5ms.

ID (mm)	Film Thickness (µm)	Length (m)	Temperature Limits (°C)	Part No.
0.1	0.1	10	-40 to 360/370	054099
0.15	1.2	10	-40 to 360/370	054106
0.15	0.25	12	-40 to 360/370	054103
0.15	0.4	12	-40 to 360/370	054107
0.15	0.25	25	-40 to 360/370	054104
0.15	0.4	25	-40 to 360/370	054108

BPX5


GC Columns and Applications

Expert Tip :

If you're having problems with solvent focusing, or early eluting peaks seem broad or lopsided in splitless injection, then try using a column with a thicker film.



ID (mm)	Film Thickness (µm)	Length (m)	Temperature Limits (°C)	Part No.
0.15	0.15	30	-40 to 360/370	054110
0.15	0.25	50	-40 to 360/370	054105
0.18	0.18	40	-40 to 360/370	054229
0.22	0.25	12	-40 to 360/370	054112
0.22	0.25	25	-40 to 360/370	054113
0.22	1	25	-40 to 360/370	054116
0.22	0.25	30	-40 to 360/370	054142
0.22	0.25	50	-40 to 360/370	054114
0.22	1	50	-40 to 360/370	054117
0.25	0.25	7	-40 to 360/370	054149
0.25	0.1	15	-40 to 360/370	0542170
0.25	0.25	15	-40 to 360/370	054100
0.25	1	15	-40 to 360/370	054121
0.25	0.1	30	-40 to 360/370	0541011
0.25	0.25	30	-40 to 360/370	054101
0.25	0.5	30	-40 to 360/370	0541025
0.25	1	30	-40 to 360/370	054122
0.25	0.25	60	-40 to 360/370	054102
0.25	1	60	-40 to 360/370	054123
0.32	1	6	-40 to 360/370	0541261
0.32	0.25	12	-40 to 360/370	054118
0.32	0.5	12	-40 to 360/370	054124
0.32	1	12	-40 to 360/370	054127
0.32	0.25	15	-40 to 360/370	054144
0.32	1	15	-40 to 360/370	054152
0.32	0.25	25	-40 to 360/370	054119
0.32	0.5	25	-40 to 360/370	054125
0.32	1	25	-40 to 360/370	054128
0.32	3	25	-40 to 350/360	054136
0.32	0.25	30	-40 to 360/370	054145
0.32	0.5	30	-40 to 360/370	0541205
0.32	1	30	-40 to 360/370	054153
0.32	0.25	50	-40 to 360/370	054120
0.32	0.5	50	-40 to 360/370	054126
0.32	1	50	-40 to 360/370	054129
0.32	0.25	60	-40 to 360/370	054146
0.32	1	60	-40 to 360/370	054154
0.53	0.25	12	-40 to 360/370	054133
0.53	1	12	-40 to 360/370	054130
0.53	3	12	-40 to 350/360	054138
0.53	0.5	15	-40 to 360/370	0541344
0.53	1	15	-40 to 360/370	054147
0.53	1.5	15	-40 to 350/360	0541347
0.53	3	15	-40 to 350/360	054159
0.53	0.25	25	-40 to 360/370	054134
0.53	1	25	-40 to 360/370	054131
0.53	3	25	-40 to 350/360	054139
0.53	0.5	30	-40 to 360/370	0541345
0.53	1	30	-40 to 360/370	054148
0.53	1.5	30	-40 to 350/360	0541348
0.53	3	30	-40 to 350/360	054160
0.53	1	50	-40 to 360/370	054132
0.53	1	60	-40 to 360/370	054158



### What is BPX5?

BPX5 is a 5% phenyl (equivalent) / 95% methyl polysilphenylene / siloxane phase.

### Why use a 5% phenyl phase?

It is estimated that 80% of all analyses could be carried out on one column with one set of dimensions and film thickness. The 5% phenyl column is still a non-polar column but the small amount of phenyl content gives the phase a little edge in selectivity compared with the 100% methyl phases.

The non-polar phases tend to have lower bleed and are more robust compared with the more polar phases.

### What can the BPX5 column be used for?

The BPX5 is a multi-purpose column which can be used for a range of applications from non-polar components (such as those in gasoline) to polar components such as amines, pesticides, triglycerides, drugs and polycyclic aromatic hydrocarbons.

For further information and details of conditions of the following chromatograms please request Publication No: AP-0069-C.

Figure 1. EPA Phenols (5ng)








Figure 2. Triglycerides in Milk Fat



**Figure 3. Analysis of EPA method 608 and additional Pesticides**



**BPX5 CAPILLARY COLUMNS ORDERING INFORMATION**

Description	Part No.	Description	Part No.
10m x 0.1mm ID - BPX5, 0.1 $\mu$ m	054099	60m x 0.22mm ID - BPX5, 0.25 $\mu$ m	054143
30m x 0.25mm ID - BPX5, 0.25 $\mu$ m	054101	15m x 0.32mm ID - BPX5, 0.25 $\mu$ m	054144
60m x 0.25mm ID - BPX5, 0.25 $\mu$ m	054102	30m x 0.32mm ID - BPX5, 0.25 $\mu$ m	054145
30m x 0.25mm ID - BPX5, 0.5 $\mu$ m	0541025	60m x 0.32mm ID - BPX5, 0.25 $\mu$ m	054146
50m x 0.15mm ID - BPX5, 0.25 $\mu$ m	054105	15m x 0.53mm ID - BPX5, 1 $\mu$ m	054147
30m x 0.32mm ID - BPX5, 0.5 $\mu$ m	0541205	30m x 0.53mm ID - BPX5, 1 $\mu$ m	054148
15m x 0.25mm ID - BPX5, 1 $\mu$ m	054121	15m x 0.32mm ID - BPX5, 1 $\mu$ m	054152
30m x 0.25mm ID - BPX5, 1 $\mu$ m	054122	30m x 0.32mm ID - BPX5, 1 $\mu$ m	054153
60m x 0.25mm ID - BPX5, 1 $\mu$ m	054123	60m x 0.32mm ID - BPX5, 1 $\mu$ m	054154
30m x 0.53mm ID - BPX5, 0.5 $\mu$ m	0541345	60m x 0.53mm ID - BPX5, 1 $\mu$ m	054158
15m X 0.53mm ID - BPX5, 1.5 $\mu$ m	0541347	15m x 0.53mm ID - BPX5, 3 $\mu$ m	054159
30m x 0.53mm ID - BPX5, 1.5 $\mu$ m	0541348	30m x 0.53mm ID - BPX5, 3 $\mu$ m	054160
15m x 0.22mm ID - BPX5, 0.25 $\mu$ m	054141	60m x 0.53mm ID - BPX5, 3 $\mu$ m	054164
30m x 0.22mm ID - BPX5, 0.25 $\mu$ m	054142		

For the complete list of BPX5 columns request AP-0069-C