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Contents

Abbreviated Instructions (detachable)

- 1 General**
 - 1.1 Important Information
 - 1.2 Safety Instructions
 - 1.2.1 Working with Rotary Vane Vacuum Pumps; Safety Information
 - 1.3 Application
 - 1.4 Design
- 2 Transport**
- 3 Installation**
 - 3.1 Base and Location
 - 3.2 Vacuum Side Connection
 - 3.3 Exhaust Side Connection
 - 3.4 Electrical Connections
 - 3.5 Monitoring Connection
- 4 Operations**
 - 4.1 Operating Medium
 - 4.1.1 Operating Medium P3
 - 4.1.2 Operating Medium F5
 - 4.2 Operating Medium Filling and Checking
 - 4.3 Starting the Pump
 - 4.4 Gas Ballast Valve
 - 4.5 Pumping off Dry Gases
 - 4.6 Pumping off Vapours
 - 4.7 Switching off the Pump
 - 4.7.1 Mode of Operation of the High Vacuum Safety Valve
- 5 Maintenance**
 - 5.1 Changing the Operating Medium
 - 5.1.1 Pump Chamber
 - 5.1.2 Motor Chamber
 - 5.2 Shutting down the Pump for Extended Periods
- 6 Technical Data**
 - 6.1 Dimension Diagram
- 7 Service Information**
 - 7.1 Customer Training
 - 7.2 Repairs in the Service Center
- 8 Trouble Shooting**

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
**Abbreviated
Instructions**
for
Rotary Vane Vacuum Pump
UNO 030B

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- 1 Filler plug, motor chamber
- 3 Vacuum connection
- 4 Exhaust connection
- 6 Filler plug, pumping system chamber
- 7 Sight glass, pumping system chamber
- 8 Operating medium drain screw, pumping system chamber
- 9 Operating medium drain screw, motor chamber (bottom)
- 10 Sight glass, motor chamber
- 13 Gas ballast valve

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These abbreviated instructions are only valid in conjunction with the full operating instruction manual.

- Please remove this sheet and place in a safe place in the vicinity of the pump.

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Initial Operation

Important:

Fill the pump with operating medium before first time operations.

Installation

- For location and surface requirements please refer to Sections 2.1 to 2.5 in the operating instructions.

Electrical Connections

- The information on the rating plate must coincide with the mains connection values.

Direction of Rotation for Three Phase Current

- Switch on briefly (max. 5 sec.).
- Place your hand on vacuum connection 3, the pump must suck.
- If the pump produces excess pressure, interchange two phases.

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Operating Medium

- Fill the pumping system chamber with the prescribed amount of operating medium P3 (supplied). The motor chamber is **already** filled with operating medium P3. **No** operating medium is supplied with pumps which are not operated with P3. In this case, the motor chamber is **not already filled**.
- Filling level: Operating medium must be visible in the upper section of the sight glass 7 when the pump is warm.

Important:

Make sure the plastic cap on the exhaust side has been removed before proceeding.

Subsequent Operation

Switching On

- The pump can be switched on in any pressure range. Maximum exhaust pressure 1.5 bar (absolute).
- The lowest run-up temperature in accordance with German Industrial Standard DIN 28 426 is + 12°C.

Checking the Operating Medium

- Check the operating medium daily when the pump is running (with vacuum connection 3 and gas ballast valve 13 closed).
- Check medium level in the upper half of the sightglass.

Gas Ballast Valve/Pumping off Vapours

- Pump off vapours **only** when the pump is warm after having been operated and the gas ballast valve is open.
- Before switching off allow the pump to go on running for approx. 30 minutes **with** open gas ballast valve so that the vapours remaining in the pump can escape and the operating medium can regenerate.

Switching Off

- Pumps can be switched off in any pressure range.
- When switching off, the integrated high vacuum safety valve blocks off the vacuum equipment from the ingress of air and also ventilates the pump.

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1 General

1.1 Important Information

Please check immediately after *unpacking* that the contents of the consignment conforms to the information given on the delivery note.

Please read the *operating instructions* before you operate the unit and follow them in all respects.

No liability will be accepted for personal injury nor material damages in the event that damage or breakdowns occur as a result of failure to comply with these operating instructions.

All units comply with the Federal German Law concerning Technical Implementations dated 24 June, 1968, Section 3. The operating instructions comply with the German Industrial Standard DIN 8418.

If you perform *repair or maintenance work* on units which have come into contact with substances which are detrimental to health, please observe the relevant regulations.

If you return units to us for repair or maintenance work, please follow the instructions below:

- *contaminated units* (radioactively or chemically etc.) must be decontaminated in accordance with the radiation protection regulations before they are returned.
- Units returned for repair or maintenance must bear a clearly visible note "*Free from harmful substances*". This note must also be provided on the delivery note and accompanying letter.
- Please use the attached attestation declaration.
- You can have the units decontaminated by us (*excepted are units with radioactive contamination*). In this case, the repair order must be extended accordingly, and the process gases with which the unit has come into contact must be stated. If this information is missing, it will be determined by us at extra cost. Any special haulage regulations must be observed.

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- We will carry out the decontamination and invoice this work to you if you have not attached the note "Free from harmful substances" to the unit or in the accompanying papers.
- "Harmful substances" are defined in European Community Countries as: "materials and preparations in accordance with the EEC Specification dated 18 September 1979, Article 2"

and in the U.S.A. as

"materials in accordance with the Code of Federal Regulations (CFR) 49 Part 173.240 Definition and Preparation".

Technical modifications reserved.

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1.2 Safety Instructions

Rotary vane vacuum pumps represent state-of-the-art technology and are designed to provide optimum operational reliability. The user may, however, be exposed to hazards if pumps are used improperly or for other than their intended purpose.

If pumps are used for any other than their intended purpose, all liability and warranty claims will lapse!

Use for the intended purpose also means that the installation, commissioning, operating and maintenance instructions of the manufacturer must be complied with.

- Operating and service personnel must have read and understood the operating instructions for the product.
- All warnings must be observed and all precautions taken.
- The operating and service personnel must be informed of all hazards which might occur in connection with the rotary vane vacuum pump. Personnel must be able to recognize dangers and take preventive measures.
- Any use by unauthorized personnel or careless handling may increase the potential danger.
- The switch-off procedures described in the attached operating instructions must be observed in all installation, commissioning, operating and maintenance work.
- No operating modes must be used which may affect the safety of the operator and the pump.
- All unauthorized modifications and alterations affecting safety are prohibited and would cause any liability or guarantee claims to be repudiated.

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- All safety protection measures must be tested on completion of electrical installation and electrical maintenance work (e.g. earthing resistance test).
- All relevant local safety and accident prevention regulations apply for operation of the system.
- Any unclear points with regard to safety, operation and maintenance should be clarified with your nearest Balzers agency or subsidiary company.

1.2.1 Working with Rotary Vane Vacuum Pumps; Safety Information

- Fill up with the prescribed amount of operating medium before initial operations.
- Never operate the pump with open vacuum flange and without strainer; particles could fall in or be sucked in and this can result in damage to the pump.
- If the pump is fixed down, a flexible flange of stainless steel or synthetic material (depending on the process) so that the rotary vane vacuum pump vibrations are not transmitted to the system or to the high vacuum pump.
- If there is no extractor on the exhaust side, the use of an oil separator is strongly recommended. The oil vapours which are generated can cause environmental pollution and would represent an infringement of the clean air laws in Germany. Balzers supplies oil separators which ensure compliance with clean air legislation.
- Only detach connection lines when the pumping stations is at a standstill so as to prevent contact with process gases.
- The pump must be fixed when it forms part of a pumping system.
- The operator is required to provide information concerning any hazards which can arise as a result of contact with dangerous media and is also required to prescribe necessary preventive measures.

1.3 Design

The pumping system and the drive motor of the Oil Immersed Single Stage Rotary Vane Vacuum Pump UNO 030 B are enclosed in a common partitioned housing. System tightness in itself is therefore of a very high level. A partitioning element, positioned between pumping system and motor, prevents dirt particles being conveyed from the pumping system side to the slide bearing motor which drives the pump. A fan prevents the pump from overheating.

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The pumping system comprises the pumping housing and the eccentrically mounted rotors with vanes within. Rotors run in slide bearings and a ball joint coupling provides the connection between pumping system and drive.

The high vacuum safety valve is controlled by a solenoid valve and seals off the vacuum chamber when the pump is not running. Simultaneously the pump is vented on the intake side.

The pump is fitted with small flange connections on the intake and the pressure side. The standard vertically mounted flanges can be converted for horizontal connection by loosening two screws.

1.4 Application

Rotary vane vacuum pumps compress gases in one operational sequence from intake pressure to atmospheric pressure. The Single Stage UNO 030 B Pump is suitable for an ultimate pressure of 10^{-2} (partial pressure).

A high level of water vapour compatibility is attained by utilising the heat loss of the inbuilt motor. Owing to its very high compression ratio this basic pump of vacuum technology encompasses a wide range of applications. It can be used in industry and research for the generation of rough and medium vacuum and also as a backing pump for vacuum pumps which do not compress against atmosphere.

2 Transport

When filled with pump fluid, pumps can be carried around within the work place by hand, using the retractable handles provided.

Caution:

Do not tilt pumps filled with pump fluid to avoid the fluid escaping through the exhaust valve.

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- 1 Filler plug, motor chamber
- 2 Carrying handle
- 3 Vacuum connection
- 4 Exhaust connection
- 5 Screw
- 6 Filler plug, pumping system chamber
- 7 Sight glass, pumping system chamber
- 8 Operating medium drain screw, pumping system chamber
- 9 Operating medium drain screw, motor chamber (bottom)
- 10 Sight glass, motor chamber
- 13 Gas ballast valve
- 19 Casing
- 120 Drive part
- 161 Terminal box


Fig. 2

Vakuumschluß 3 horizontal
Vacuum connection 3, horizontal
Raccord vide 3, horizontale

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3 Installation

Warning!

Fill the pump with operating medium prior to putting it into operation for the first time (see Section 4.1).

3.1 Base and Location

Fig. 1

- Place the pump on an even, flat surface.
- It is unnecessary to fix the pump to the surface.
- Sight glasses 7 and 8 must be visible.
- Maximum ambient temperature + 40 °C.
- Ensure adequate air circulation if installing in a closed area. The intake opening of the fan must be kept free.

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3.2 Vacuum Side Connection

Fig. 1, Fig. 2

- The connection between recipient and pump should be as short as possible.
- Vacuum connection 3 is mounted vertically and can be turned horizontally (Fig. 2) if required.
- Any angle of connection is possible by using intermediate flanges (see Section 11, Accessories).
- Welded piping must be descaled and any foreign particles must be removed before installation.

3.3 Exhaust Side Connection

Fig. 1

- Lay exhaust lines without built-in shut-off devices.
- If a shut-off valve is installed for internal operating reasons or if there is a risk of excess pressure developing in the lines, official accident prevention regulations must be complied with.
- The exhaust line diameter should not be less than 25 DN.
- The line from the pump must be laid with a downward gradient so that condensate is prevented from flowing back into the pump. Otherwise a trap must be installed.

Important:

Expelled gases and vapours can represent a hazard to health and/or the environment.

For this reason the relevant regulations must always be complied with when working with toxic substances. Only officially approved filters may be used for the separation of such substances. This also applies when handling toxic or chemically unstable gases and to gases susceptible to polymerisation or to the formation of peroxide (see Section 11, Accessories).

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Fig. 4 to Fig. 7

Please note:

Electrical connections are subject to whatever local regulations are in force.

- Mains power voltage and frequency must comply with the voltage and frequency data shown on the rating plate.

Important:

- Check the direction of rotation of the pump before filling up with operating medium.
- Remove plastic sealing cap on the exhaust side.

For three phase drives, check the direction of rotation as follows:

- Switch on briefly.
- Place your hand on vacuum connection 3, the pump must suck.
- If the pump produces excess pressure, interchange two phases of the connecting cable.

Pumps with an **alternating current drive** are always supplied with an ON/OFF button, unless otherwise specifically ordered. To monitor the temperature, a switch contact is incorporated in each motor winding, with which a contactor suitable for the motor in question can be switched.

Pumps with **three phase current drive** can be supplied with or without a motor protection switch, as stipulated in the order. To monitor the temperature, three PTC resistors are incorporated in each motor winding. To switch the motor contactor, a suitable switching unit is needed when connection is executed on site.

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Fig. 4, Fig. 7

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3.5 Monitoring Connection

Fig. 3, Fig. 12

- The pumps are equipped for the connection of a Rotation Speed/Operating Hours Control Unit RCI 002. Installation and use of this unit are described in Operating Instructions PK 800 074 BD,E,F.

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- The unit is connected to the terminals ~~/ on the motor terminal board.
- Feed-throughs are provided in terminal box 161 for cable screw connections Pg. 7.

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Alternating current drive
Ready for connection, with ON/OFF switch at the terminal box
cover temperature monitor, cable and plug.

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Alternating current drive
ON/OFF switch provided by customer.

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Three-phase. Customer connection, with contactor control and tripping unit TMA 4R for temperature monitoring.

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Three-phase, with ON/OFF switch at the terminal box cover and excess current tripper.

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- 12 Motor contactor with delayed cut-out
- 131 Rectifier PCB
- 141 Cover, terminal box
- 145 Mains connection
- 150 Terminal board
- 161 Terminal box
- 168 Excess temperature cut-out

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4.1 Operating Medium

4.1.1 Operating Medium P3

Standard model pumps are filled with operating medium P3 for testing purposes. P3 oil is distinguished by its low vapour pressure, good lubricating action, favourable alkaline reserve and high oxidation resistivity.

No guarantees regarding the attainment of ultimate pressures, and oil-related damages apply if oil other than P3 is used.

Operating medium P3 must be disposed of in accordance with local regulations.

4.1.2 Operating Medium F5

Only the synthetic operating medium F5 may be used for corrosive gas model pumps. The use of other operating media must be cleared by the manufacturer.

Caution:

When working with synthetic operating media, the manufacturer's directions for use must be observed. If F5 is heated to over 250°C, poisonous gases damaging to the respiratory system are produced. F5 must not be allowed to come into contact with tobacco products. When handling chemicals, the relevant manufacturer's instructions must be followed.

Operating media which have been enriched with toxic or corrosive substances must be specially handled and all pertinent regulations must be observed.

4.2 Filling and Checking

Fig. 1

Pump chamber

- Unscrew filler screw 6.
- Pour in the operating medium supplied. Quantity: see rating plate.
- Screw in filler screw 6; pay attention to the O-ring!
- Check the operating medium level only with the pump running (warm).
- Close vacuum connection 3 and gas ballast valve 13.
- Filling level: upper half of sight glass 7. Minimum level: middle of sight glass 7.
- The operating medium can be topped up with the pump running at ultimate vacuum level.
- For continuous operation, check daily, otherwise, check every time the pump is switched on.

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Motor Chamber

Note:

On standard model pumps the motor chamber has been filled with operating medium P3 in the factory.

- Unscrew filler screw 1.
- Pour in operating medium until the level reaches the middle (the maximum) of sight glass 10.
- Screw in filler screw 1 keeping an eye on the O-ring.
- A periodic check is sufficient.
- When the motor is warm, the operating medium can ascend to the upper half of sight glass 10.

4.3 Starting the Pump

- The pump can be started in any pressure range.
- Minimum start-up temperature, in compliance with German Industrial Standard DIN 28 426 + 12 °C. Other temperatures are subject to agreement.

4.4 Gas Ballast Valve

Fig. 8

Gas ballast valve 13 is closed in position 0 and open in position I. An intermediate setting is not possible.

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- 3 Vacuum connection
- 4 Exhaust connection
- 13 Gas ballast
- 14 Dosing screw
- 40 Cover with position indicator

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4.5 Pumping off Dry Gases

No special measures are necessary. The gas ballast valve can remain closed.

4.6 Pumping off Vapours

To prevent their condensation in the pump, vapours should only be pumped off with the pump at operating temperature and the gas ballast valve 13 open. Observe pump water vapour tolerance levels (see Section 1, Technical Data). When the process has been completed, allow the pump to continue running for approx. 30 minutes for operating medium cleaning purposes.

4.7 Switching off the Pump

The pump can be switched off in any pressure range. On switching off, the built-in high vacuum safety valve 21 to the vacuum apparatus closes automatically and vents the pump.

Note:

To switch off the alternating current model pump, the button must be held depressed approx. 2 seconds.

4.7.1 Mode of Operation of the High Vacuum Safety Valve Fig. 9

The automatic high vacuum safety valve 21 is open when the apparatus is vented. When the pump is switched on, voltage is generated over the coil 99 with permanent magnet 100, solenoid 174 is activated and attracts solenoid core 173. Air intake 17 is then closed and plunger 26 remains in the open position due to the pressure of spring 33.

If the rotor 108 stops rotating when the equipment is in its evacuated state (power failure or pump switched off), the coil with the permanent magnet no longer transmits any voltage to solenoid coil 174 and the spring pressure opens solenoid core 173. The air which now flows in, presses valve plate 34 against the vacuum side 3 by way of valve plunger 26 and closes off the equipment vacuum-tight. At the same time, pumping system 194 is vented.

After the vacuum pump has been switched back on again, the air intake is shut off by solenoid core 173 and the interior of the valve is evacuated. When the pressure in the valve has decreased enough so that the pressure of spring 33 is sufficient to push valve plunger 26 against the pressure prevailing in the equipment, the valve opens.

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6. Maintenance

When carrying out servicing work, take the following precautions:

Always ensure the pump cannot be switched on when carrying out any work on the pump. If necessary, remove the pump from the system to carry out inspection work.

- ➔ Only dismantle the pump as far as is necessary in order to repair defects.
- ➔ Dispose of used operating fluid in compliance with local regulations.
- ➔ When working with synthetic operating media, toxic substances and substances contaminated with corrosive gases, the relevant instructions governing their use must be observed.
- ➔ Only use benzene or similar agents for cleaning pump parts. Do not use soluble cleaning agents.

5.1. Changing The Operating Fluid

The rate of deterioration of the operating fluid depends on the pump applications.

- ➔ The level of deterioration of organic operating fluids (for example P3) can be read off from the colour scale in accordance with DIN 51578 on the supplementary sheet PK 800 219 BN/B, on request.
 - The colour scale enables precise colour determination.
- ➔ Fill a specimen in a test tube or some similar vessel and test by holding against the light.
- ➔ Where the discolouration is dark yellow to red brown (equivalent to 4 ... 5 on the scale) change operating fluid.

Where organic operating fluids, for example P3, are involved, the operating fluid should be changed at least once a year!

The operating fluid temperature can be as high as 80 °C. During maintenance and repair work, process related toxic gases and vapours can escape from the operating fluid which may become contaminated with harmful substances (radioactive, chemical etc.).

Disposal of used operating fluid is subject to the relevant regulations.

Operating fluid can contain substances from the medium pumped. Operating fluid must be disposed of in accordance with the local respective regulations.

Please request safety instruction data sheets for operating fluids from PFEIFFER.

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5.1.1 Pump Chamber

- Switch off the pump.
- Unscrew drain screw 8 and drain off the operating medium.

Caution:

The operating medium temperature can be as high as 100°C. During maintenance and repair work application related toxic gases and vapours can escape from the operating medium which may become enriched with harmful substances (radioactive, chemical etc.).

When working with F5 it is imperative to observe the precautionary measures referred to under Section 4.1.2. The disposal of used operating medium is subject to the relevant regulations.

- Screw in drain screw 8; pay attention to the position of the O-ring.
- Let the pump run for approx. 10 seconds with the vacuum connection open.
- Drain off the remaining operating medium.
- If very dirty, change the medium several times.
- Fill with operating medium and check the level of medium as described in Section 4.2.

5.1.2 Motor chamber

Change the operating medium every 10,000 operating hours.

- Unscrew drain screw 9 under drive component 120.
- To open to air, unscrew filler screw 1.
- Drain off the operating medium.
- Screw in drain screw 9; pay attention to the O-ring.
- Fill with operating medium and check as described in Section 4.2.

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5.2 Shutting down the Pump for Extended Periods

Fig. 1

If the pump is shut down over an extended period of time, the entire pumping system must be adequately protected against corrosion.

- Switch off pump.
- Remove operating medium drain screw 8 and drain operating medium as described in Section 5.1.

Note the important advice given in Section 5.1!

- Insert drain screw 8; check the O-ring seating.
- Allow the pump to run for approx. 10 seconds with open vacuum connection.
- Drain remaining operating medium.
- Fill up with 3.5 l new operating medium.

If the operating medium quantities shown in the table are used, the liquid level is extremely high, but this is necessary to protect all pump parts against corrosion.

Caution:

Before putting the pump back into operation, the operating medium level must be reduced to normal.

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7 Service information

6.1 Customer Training

We hold regular seminars on vacuum technology in our works. These seminars contain a theoretical and a practical part. The theoretical part deals with the design and function of the respective pump. In the practical part, pumps are dismantled and explanations given on how defective pumps are repaired. If you are interested in participating in a course we would be pleased to send you further information.

Should you have any questions, please contact our sales companies or agencies which will be pleased to assist you.

Please see back cover for addresses.

6.2 Returning to Service Center for Repair

We request that you contact your local Pfeiffer Service Center should problems arise.

All repair orders are carried out on the basis of our General Terms of Delivery only. For repairs, the warranty period is 6 months. Repair work is performed without checkback if the repair costs do not exceed 50 % of the nominal value of the article or DM 5,000.00, whichever is less.

*Please note:
The operating medium in the pump must be drained before you return the pump for repair.*

When returning for repairs, the important introductory notes contained in these instructions must be observed.

Recommendation:
We recommend that you take out a service contract with your local Pfeiffer Service Center, thereby avoiding possible extended periods of inactivity which could be caused by major damage to the pump.

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8 Trouble Shooting


Trouble shooting is divided into five main categories, depending on the source of malfunction in the pump. These sources of malfunction can be ascribed to a number of possible causes. The measures necessary to correct the malfunction depend on the cause.


The references to 9.1, 9.2 etc. show where the instructions concerning assembly and dismantling can be found.

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



Pump starts but stops again after a short time

Is the fan also turning

no → Check fan and replace if necessary

yes →

Three phase current. Is the overload circuit breaker set correctly

nein → Check setting

yes →

Inspection only to be performed by qualified personnel

yes →

Is the pumping system clean

no → Clean pump system, see Section 9.11

yes →

Is the coupling in order

no → Check coupling, see Section 9.11

yes →

Is the slide bearing in the pumping system in order

no → Check slide bearing, see Section 9.12

yes →


Is the slide bearing in the motor in order

no → Check slide bearing, see Section 9.12

yes →

Inform **PFEIFFER**-Service, see Section 7.2

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9.1.1 Alternating Current Model

Wiring diagrams PK 212 824 -S and
PK 212 825 -S.

- Open cover 141 of terminal box 161.
- Switch on the pump.
- Measure the voltage at L1-N and U1-U2. If voltage is present at L1-N but not at U1-U2, ON/OFF switching is faulty.
- Check fuse F2 on the rectifier PCB 131.
- If F2 is o.k., either the motor contactor or the temperature monitoring is faulty.
- Check the temperature monitoring as described in Section 9.3.
- If the temperature monitoring is o.k., contactor 12 must be exchanged complete with the delayed cut-out. The delayed cut-out bridges any interruption in the mains supply for up to max. 1 second.

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Fig. 11, Fig. 12

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
Legend, Wiring Diagrams
PK 212 824 -S and PK 212 825 -S

- 11 LED display
- 12 Contactor with delayed cut-out
- 131 Rectifier PCB
- 132 Run-up control
- 141 Terminal box cover
- 145 Mains connection, pump
- 146 Mains connection, RCI 002
- 150 Terminal board
- 161 Terminal box
- 162 Rotation speed monitor RCI 002
- 166 Frequency change
- 168 Excess temperature cut-out

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
Alternating current motor with ON/OFF
 switching on site

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Alternating current motor, ready for connection,
 with ON/OFF switching with terminal box
 cover, temperature monitoring, cable and plug

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9.1.2 Three-Phase Model

Wiring diagrams PK 212 826 -S
and PK 212 827 -S

- Open cover 141 of terminal box 161.
- Switch on the pump.
- Check the voltage of L1, L2 and L3 at the fuses 10A to U1, V1, W1 at the motor terminal board. Exchange any defective parts.
- Check whether voltage is present at L and at terminal 14 of the switching unit 170. If not: Check temperature monitoring as described in Section 9.3.

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
Legend for Wiring Diagrams
PK 212 826 -S and PK 212 827 -S

- 11 LED display
- 12 Contactor with delayed cut-out
- 131 Rectifier PCB
- 132 Run-up control
- 141 Terminal box cover
- 145 Mains connection, pump
- 146 Mains connection, RCI 002
- 150 Terminal board
- 161 Terminal box
- 162 Rotation speed monitor RCI 002
- 169 TF-PTC resistor
- 170 Tripping unit for 169

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
Three phase motor with ON/OFF switching
 in terminal box cover and overload circuit
 breaker

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Three phase motor for connection on site
 with contactor control and tripping unit for
 temperature monitoring

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9.2 Voltages and Frequencies

9.2.1 Alternating Current Motor

The motor is designed and switched in accordance with the specifications given in the order. Voltage and frequency are given on the rating plate.

It is possible to change the connection of the motor in accordance with the required voltage and to mark the rating plate accordingly.

Two motors are used for the alternating current drive:

- A – for 230 V ± 15 % / 50 Hz
and 208 V ± 10 % / 60 Hz.
- B – for 100 V ± 10 % / 50 Hz
and 110 V ± 15 % / 60 Hz.

Both motors are designed for 50 and 60 Hz. Wiring for the respective voltage in accordance with the table. See Section 9.5 for switching the run-up control frequency.

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9.2.2 Three-Phase Motor

The three-phase motors used are designed for 50 and 60 Hz. A frequency switch-over is not required. Motor is wired in accordance with the table.

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9.3.1 Alternating current Motor

Fig. 13

- The temperature-dependent switch contact in the motor winding is closed at temperatures below 130°C and opens when the temperature exceeds 130°C, terminal 20/21.
- If the switch is faulty, the drive component 120 must be exchanged (Fig. 12).

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Switching: alternating current motor with excess cut-out 20/21

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Switching: three-phase motor with tripping unit TMA 4R for excess temperature cut-out via three PTC resistors in the motor winding.

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9.3.2 Three-Phase Motor

Fig. 11, Fig. 14

The switch contact in tripping unit 170 is closed when the pump is cold and opens at a pump temperature of approx. 130°C.

Checking:

- Switch off the pump and let it cool down.
- Disconnect the cable from terminals 10 and 11 on the motor terminal plate 150.
- Measure the resistance at terminals 10 and 11. At approx. 20 to 30°C this should be 200 – 250 Ω. If the resistance is not measurable, drive component 120 must be exchanged.
- If the PTC resistors are in working order, measure tripping unit 170. For measurement purposes, the PTC resistors must be connected and voltage must be present at the tripping unit 170. If voltage is present at terminal 13 but not at terminal 14, the tripping unit 170 must be exchanged.
- Check the operating rotation speed. Measure at the terminal ~/~: At a mains frequency of 50 Hz the frequency should be approx. 190 Hz and at a mains frequency of 60 Hz, it should be approx. 230 Hz. If the frequency is lower, the pumping system or the bearing in the drive component is sluggish. Dismantle and clean the pump, exchange any faulty parts as described in Section 9.11 and 9.12.

9.4 Motor Winding

9.4.1 Alternating Current Motor

Fig. 12, Fig. 13, Fig. 15

- Disconnect the pump from the mains.
- Unscrew terminal box cover 141 and open it.
- Unscrew all U1, U2, U7, U8, U9 and Z2 connections.
- Check the motor windings as shown in Fig. 15 for winding shorts and shorts-to-earth.

If there is no passage through any of the motor windings or if there is a short-to-earth, the drive component 120 must be exchanged.

9.4.2 Three-Phase Motor

- Disconnect the pump from the mains.
- Unscrew terminal box cover 141 and open it.
- Unscrew all U1, V1, W1 and U2, V2, W2 connections.
- Check the motor windings for shorts and earthing.

If there is no passage through any of the motor windings or if there is a short-to-earth, the drive component 120 must be exchanged.

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Fig. 11, Fig. 14

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Location of motor windings on the terminal board of the alternating current motors.

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- Switch off the pump.
- Unscrew the run-up control 132 and lean it against the pump for inspection.
- Measure the passage through spark quenching resistor R 505.
- Switch the pump on.

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- Double the mains voltage must always be present at the operating capacitor CB 135. The starting capacitor CA 134 only has double the mains voltage during the run-up phase. After the switching point has been reached (at a mains frequency of 50 Hz, the switching frequency is 176 to 180 Hz and at a mains frequency of 60 Hz, the switching frequency is 211 to 216 Hz), the mains voltage drops to zero. If the voltage deviates, the respective capacitor must be exchanged.
- Measure the input voltage at terminal ~/~ 4 - 4.5 V alternating current must be present.
- If you cannot find the defect, the run-up control 132 must be exchanged.

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

Frequency switching from 50 to 60 Hz or vice-versa is performed on the PCB using the toggle switch unit 166 at switch 1 (position "OFF" 50 Hz).

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- Check operating rotation speed. Measure at terminal ~/~: At a mains frequency of 50 Hz, the frequency is approx. 190 Hz and at a mains frequency of 60 Hz, approx. 230 Hz. If the frequency is lower, the pump system or the motor bearing is sluggish. Dismantle and clean the pump, exchange any faulty parts as described in Section 9.11 and 9.12.
- If the starting capacitor CA 134 does not switch off with a one-phase motor, the switching point for the starting capacitor CA 134 must be checked between the test point TP and the earthed conductor 196. At a mains frequency of 50 Hz, the switching frequency is 176 to 180 Hz and at a mains frequency of 60 Hz, the switching frequency is 211 to 216 Hz. If necessary, readjust the switching point at the switches 2 to 4 of toggle switch unit 166.

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9.6 Resetting the Silencer

Fig. 18

If the silencer is not functioning, you will hear a knocking noise at the end vacuum (oil thumping). This does not cause damage to the pump but it may be disturbing in quiet rooms. Setting is executed with the pump running, the intake socket closed and the gas ballast valve closed, using a normal screwdriver. The pump must have reached its operating temperature.

- Remove cover 40 from gas ballast valve 13.
- With a screwdriver, open dosing screw 43 approx. 2 turns (flushing) and then fully close it again (clockwise).
- Wait until the knocking noise of the oil is audible.
- Open dosing screw 43 slowly (1/4 to 1/2 a turn, depending on the type of pump) until the knocking noise is no longer audible.
- Refit cover 40 onto gas ballast valve 13.

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Fig. 18

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3 Vacuum connection
4 Exhaust connection
13 Gas ballast valve
40 Cover with position indicator
43 Dosing screw

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9.7 Gas Ballast Valve

Fig. 19

Gas ballast valve 13 can only get dirty if the intaken ambient air is dusty.

- Unscrew screws 39 and pull gas ballast valve 13 out of the pump housing.
- Loosen hose clip 53 and take off hose 54.
- Remove cover 40.
- Take off sleeve 42, paying attention to O-ring 41.
- Screw off dosing screw 43.
- Screw off connecting piece 48, paying attention to the O-ring 45.
- Unscrew screw 50 and remove valve spring 47.
- Clean all parts.
- Reassemble the gas ballast valve in reverse order. Valve spring 47 must be screwed to part 46, as shown in detail "A" in Fig. 20.
- Plug hose 54 onto part 48, paying attention to the series of compression springs 147 in hose 54.
- Fasten hose clip 53.
- Fit gas ballast valve 13 in pump housing 120. Make sure that O-ring 121 is in the correct position.

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Fig. 19

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Detail "A"
- Insertion position of valve spring 47

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- 13 Gas ballast valve
- 39 Screw
- 40 Cover with position indicator
- 41 O-ring
- 42 Sleeve
- 43 Dosing screw
- 44 O-ring
- 45 O-ring
- 46 Valve housing
- 47 Valve spring
- 48 Connecting piece (hose)
- 49 Spring sleeve (dowel pin)
- 50 Screw
- 51 Gas inlet
- 53 Hose clip
- 54 Hose
- 120 Drive part

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9.8 Dismantling and Assembling the Casing

Fig. 20

- Drain off the operating medium as described in Section 5.1. If only the pump chamber is to be inspected, it is not necessary to drain off the operating medium from the motor chamber.
- Unscrew nuts 107 and remove washers 167.
- Take off casing 19 in axial direction, paying attention to O-ring 18.
- Clean the interior of casing 19.

Casing 19 is reassembled in reverse order. Ensure that the O-ring 18 is in the correct position.

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Fig. 20

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9.9 High Vacuum Safety Valve

Fig. 21, Fig. 22

- Switch off the pump.
- As the high vacuum safety valve 21 shuts off at the vacuum side, the vacuum side of the pump must also be vented.
- Drain off the operating medium as described in Section 5.1.
- Dismantle casing 19 as described in Section 9.8.
- Disconnect connector 20.
- Loosen retainer 15.
- Pull off high vacuum safety valve 21 and dismantle it as shown in Fig. 22.
- Do not dismantle parts 26 – 28 and 31 – 34 for cleaning purposes.
- Clean the parts.
- Check the parts for wear; exchange any defective parts.
- Reassemble the parts in reverse order.
- Insert angle plate 36 (Fig. 26).
- Insert O-ring 122. Angle plate 36 must protrude approx. 2 mm above the O-ring 122.
- Slightly tighten clamping screw 37, to hold angle plate 36 in its present position.

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Fig. 21, Fig. 22

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- 15 Retainer
- 16 Screw
- 17 Stud bolt
- 18 O-ring
- 20 Connector
- 21 High vacuum safety valve
- 22 Solenoid valve
- 23 Screw
- 24 Cover
- 25 O-ring
- 26 Valve plunger
- 27 Valve plate
- 28 O-ring
- 30 O-ring
- 31 Circlip
- 32 Valve spindle
- 33 Pressure spring
- 34 Valve plate
- 35 Valve housing
- 129 Cable

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Fig. 23, Fig. 24

Function

In a currentless state, solenoid valve 22 is open. Activation is executed via the compression springs on the solenoid core 173. If solenoid coil 174 is supplied with current, the core is drawn to the valve seat and the valve closes (Fig. 9).

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Fig. 23, Fig. 24

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20 Cable connection
171 O-ring
173 Solenoid core with compression spring
174 Solenoid coil
175 Valve housing
176 Nut

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- 1 Filler plug, motor chamber
- 20 Connector
- 22 Solenoid valve
- 54 Hose
- 99 Coil
- 100 Permanent magnet
- 101 Pole shoe
- 118 O-ring
- 128 O-ring
- 129 Cable
- 131 Rectifier PCB
- 141 Cover, terminal box
- 147 Series of pressure springs
- 150 Terminal board
- 161 Terminal box
- 177 Cable feedthrough
- 178 Terminals, control voltage

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9.10.1 Checking the solenoid valve voltage
Fig. 25

If solenoid valve 22 is not working, the following test should be carried out before dismantling:

- Screw off cover 141 from terminal box 161.
- Switch on the pump.
- Measure the voltage at rectifier PCB 131:
Between $\sim/\sim = 4$ to 4.5 V
Between $+/- = 2.8$ to 3.5 V
- If voltage is present, the fault is located either in solenoid valve 22 or in the cables between PCB 131 and solenoid valve 22.

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- If no alternating current voltage is present, the fault can be found in coil 99 or in the permanent magnet 100.
- If no D.C. voltage is present, the fault lies in the rectifier PCB 131.

Please note:

A dirty solenoid valve 22 can have the same effect as a mains failure.

Important

All screwed parts of solenoid valve 22, including those on the valve cover 24, (Fig. 22) may not be subjected to a tightening pressure in excess of 4 Nm.

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

- 11 LED display, "Pump running"
- 22 Solenoid valve
- 99 Coil
- 100 Permanent magnet
- 101 Pole shoe
- 131 Rectifier pcb
- 150 Terminal board
- 161 Terminal box
- 178 Terminals, control voltage
- C1 Capacitor
- n1 Rectifier

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9.10.2 Dismantling and assembling the solenoid valve

Fig. 22, Fig. 23, Fig. 24

Please note:

The operating medium must be drained from the pumping system before dismantling solenoid valve 22.

- Dismantle casing 19 as described in
Section 9.8.
- Disconnect both cables 129 at connec-
tor 20.
- Unscrew the solenoid valve 22 from
cover 24.
- Take off O-ring 30.
- Pull valve housing 175 out of the sole-
noid coil 174.
- Screw nut 176 off valve housing 175
and pull the solenoid core 173 out of
part 175.
- Do not remove the compression spring
from solenoid core 173.
- Clean all parts and check them for
wear.
- Solenoid valve 22 can only be exchan-
ged as a complete unit.
- Reassemble the solenoid valve and
pump in reverse order.

9.10.3 Checking the Cable Feedthrough Fig. 24

- Dismantle casing 19 as described in
Section 9.8 and high vacuum safety
valve 21 as described in Section 9.9.
- Disconnect both cables 129 from con-
nector 20.
- Check cables 129 for passage of cur-
rent and shorting.
- If a cable is faulty, the cable feed-
through (bushing) 177 must be
exchanged as a complete unit.
- Unscrew filler screw 1 (motor cham-
ber), paying attention to O-ring 118.
- Unscrew cover 141 from the terminal
box 161.
- Disconnect cable 129 from the termi-
nal board 150 and terminal box 161.
- For alternating current models,
dismantle run-up control 132 (Fig. 12).
- Screw off terminal box 161.
- Disconnect cable 129 from the termi-
nal board 150.
- Unscrew cable feedthrough 177 and
pull it out of the drive component 120
together with cables 129.
- Fit new cable feedthrough 177 with
cables 129.
- Reassemble the remaining parts in
reverse order.

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9.11 Pumping System

9.11.1 Disconnecting the Pumping System from the Drive Component Fig. 26

The high vacuum safety valve can maintain its shut-off function when the pumping system is disconnected. It need not be dismantled.

- Dismantle casing 19 as described in Section 9.8.
- Dismantle high vacuum safety valve 21 as described in Section 9.9, or simply screw off retainer 15.
- Unscrew nut 83, paying attention to spring washer 183.
- Remove pumping system in axial direction.
- Pumping system as applicable, out of the housing until it is possible to loosen hose clip 53 and pull off hose 54.
- Take off the pumping system.

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9.11.2 Fitting the Pumping System to the Drive Component

Fig. 11, Fig. 12, Fig. 26

- Place the drive gib of coupling half 114 in a horizontal position and place on coupling disk 55.
- Place rotor carrier key in vertical position.
- Carefully push the pumping system onto the two locking screws 52.
- Connect gas ballast hose 54 and fasten it with hose clip 53. Pay attention to the series of compression springs in the hose.
- Push the pumping system against lantern 119.
- Push on spring washer 183.
- Screw on nuts 83, align the pumping system horizontally and tighten the nuts to 4 Nm.
- Lightly oil O-ring 72 and insert connecting piece 70 into the pumping system.
- Align flat connecting piece 70 horizontally.
- Lightly oil O-ring 71 and insert it.
- Check the position of O-ring 123 and oil the shaft of angle plate 36.
- Push high vacuum safety valve 21 onto angle plate 36.
- Fasten the high vacuum safety valve 21 onto connecting piece 70 with retainer 15.
- Tighten clamping screw 37.
- Fit casing 19, ensuring that O-ring 18 is in the correct position.
- Fit vacuum fitting 3, ensuring that O-ring 122 is in the correct position.
- Fill the pumping system chamber with operating medium, as described in Section 4.2.
- For pumps with three phase drives, check the direction of rotation as described in Section 3.4.
- Fit the pump into the plant.

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9.11.3 Dismantling the Pumping System
Fig. 27

Please note:

The sealing surfaces of all components of the pumping system have either been lapped or ground. Take care not to damage these surfaces when dismantling or assembling the systems.

- Clean the external surfaces of the pumping system.
- Pull connecting piece 70 out of cylinder 57.
- Dismantle exhaust valve 195 (3 pieces).
- Unscrew screw 61. Pay attention to spring washer 189.
- Dismantle the pumping system in axial direction.
- Clean all parts.

The rotor is borne eccentrically in the cylinder. In cylinder 57, there is a recess or limit stop between the intake and exhaust channel. The remaining air gap between the cylinder and the rotor is sealed by the operating medium.

- Check whether the rotor has run into this recess.
- Check vane springs 59.
- Check whether it is easy to move the vane 58 in the rotor 62.
- Check whether the rotor pivot has run into the sliding bearings of parts 56 and 60.
- Blow out the operating medium channels with compressed air.
- Check pressure valve 63. This valve can only be exchanged as a complete unit.
- Examine coupling disks 55 and coupling half 114 for wear.
- Check the coupling driver at the rotor for wear.
- Exchange any defective parts.

Parts 56, 57, 60 and 62 are not exchangeable. The whole pumping system must be exchanged. All other parts can be exchanged individually.

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9.11.4 Assembling the Pumping System

- The pumping systems are assembled in reverse order to that in which they were dismantled.
- Oil all pump components prior to assembly.
- Screws 61 and 92 must be tightened with a torque spanner to 4 Nm. Excessive tightening of these screws will distort and jam the pumping system.
- Prior to proceeding with final assembly of the individually assembled pumping system, check that the rotor rotates in the respective pumping system.
- Fill with operating medium as described in Section 4.2.

Warning!

Operating medium boreholes can be found in one of the two bearing pins of the rotor 62, Fig. 27. When assembling, ensure that this bearing pin is always in support plate 56.

9.12 Drive

Fig. 11, Fig 12, Fig. 28

Warning!

Disconnect the pump from the mains.

9.12.1 Dismantling the Drive

See Section 10, Tools for the tools required for assembly and dismantling.

- Dismantle the pump as described in Section 9.11.
- Drain off the operating medium out of the motor chamber as described in Section 5.1.2.
- Place the pump vertically on bearing cover 102 for the three phase model or on run-up control 132 for the alternating current model.
- Unscrew screws 112.
- Pull the rotor 108 with lantern 119 out of the drive component 120 at the stud bolts 52. If the lantern 119 is stuck, unscrew stud bolts 52 and push it out with cheesehead screws in place of the stud bolts.

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- Do **not** dismantle centering sleeve 81.
- Unscrew screw 113 and screw in a hexagonal socket screw M4 x 35 in its place. Then pull off coupling half 114 with the two-armed puller 1 (see Section 10).
- Dismantle parts 109, 110, 115, 116, 156 and 157.
- Set the pump housing on its feet.
- For the alternating current model, dismantle run-up control 132.
- Unscrew cover 141 from terminal box 161.
- Disconnect the cable from terminal board 150 and terminal box 161.
- Dismantle terminal box 161 and bearing cover 102.

Caution!

Coil 99 is electrically wired to terminal board 150.

- Unscrew the cable for coil 99 at terminal board 150.
- Remove coil 99 from the bearing cover 102. Clean all parts, with the exception of coil 99, in petroleum ether or similar.
- Check parts 109, 110 and the two sliding bearings in parts 102 and 119, and the sliding bearing raceways at the rotor 108 for wear.
- Exchange any defective parts.
- If the sliding bearings have run into the raceways, parts 102 and/or 119 must be exchanged.
- Only exchange parts 109, 110 and 116 all together.

Exchange parts

As component 120 according to the spare parts list.

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9.12.2 Assembling the Drive

- Fit coil 99 in bearing cover 102. Circlip 106 must fit properly into the slot.
- Connect the cable of coil 99 to terminal board 150 as described in Section 9.10.
- Screw on bearing cover 102.
- Place drive component 120 (in effect the pump housing) onto bearing cover 102.
- Screw in both screw plugs 1 and 9 with O-rings 118 and 117.
- Stand the rotor 108 upright and completely assemble it.
- Insert axial face seal 109 into the mounting tool 8 (see Section 10, Tools) and push it into lantern 119.
- Insert O-ring 111 and push the lantern onto the rotor shaft.

Assemble sliding ring 110 as follows:

- Push mounting sleeve 4 (see Section 10, Tools) onto the rotor shaft until it reaches the limit stop and lightly oil it.
- Push sliding ring 110 together with O-ring 116 over the mounting sleeve 4 onto the rotor shaft, until it reaches the limit stop.
- Fit supporting ring 156 and circlip 157.
- Fit key 115 and coupling half 114 and secure with screw 113.
- Screw stud bolts 52 into the lantern 119 and insert the rotor in the pump housing. The stud bolts must be parallel to the base.
- Screw on lantern 119 with screws 112.
- Stand drive component 120 on its feet.
- Fit terminal box 161 complete; pay attention to the profile rings 97.
- For alternating current models, fit run-up control 132.
- Connect the cables to the terminal board in accordance with the wiring diagram.

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10 Tools

Fig. 30

The assembly and mounting tools are not included in the normal consignment; they can be ordered as a boxed set under order nr. PK 002 985 -T or separately under the order numbers given.

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- 1 Two-armed puller PK 002 976 -R
(to pull coupling 114 off the rotor shaft)
- 2 Mounting clamp PK 002 974 -X
(to dismantle and assemble the pressure relief valve 190)
- 3 Face spanner PK 002 977
(to dismantle and assemble the sight glass)
- 4 Mounting sleeve PK 002 975
(to fit the sliding ring 110)
- 5 T-handle P 0984 226
- 6 Extension P 0984 225
- 7 Drift punch PK 002 979
(to dismantle the axial face seal 109, see Fig. 31)
- 8 Mounting tool PK 002 978 and handle PK 002 959
(to fit the axial face seal 109, see Fig. 32)
- 9 Spanner heads SW17
(to assemble and dismantle the cable feedthrough 177)

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9 Mounting tool with handle
10 Drift punch
109 Axial face seal
119 Lantern

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11 Accessories

A comprehensive range of accessories can be found in Section A of the Balzers catalogue, to: a) extend the application range of the pumps and b) protect the pump from wear, when used under extreme operating conditions.

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12 Spare parts

Please note!
When ordering spare parts please be sure to state the model and its code number.

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¹⁾ Axial face seal 109 must be exchanged together with sliding ring 110 and O-ring 116

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

¹⁾ Not included in the consignment

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
Wechselstrom
 Alternating current
 Courant alternatif

Drehstrom
 Three phase current
 Courant triphasé

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Safety information for the return of vacuum systems and components

General Introduction

The employer (user) is responsible for the health and safety of his employees. This also applies to all those persons who come into contact with the vacuum components either at the user's, or manufacturer's premises during repair or service. The contamination of vacuum systems and components has to be declared and the Health and Safety Declaration form completed.

Health and Safety Declaration

Those persons carrying out repair or service have to be informed of the condition of the components. This is the purpose of the 'Declaration of Contamination of Vacuum Equipment'.

This declaration must be sent directly to the manufacturer or representative company. A second copy must be sent with the equipment.

Despatch

When returning equipment the procedures set out in the Operating Instructions must be followed. For example:

- drain the pump operating medium,
- neutralise by flushing with gas,
- remove filter elements,
- seal all outlets,
- seal in heavy duty polythene or a bag,
- despatch in suitable transport container.

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⇒ **GB, IE**

Declaration of conformity pursuant to the following EU directives:

- **Machinery Directive 89/392/EEC**
- **Electromagnetic Compatibility Directive 89/336/EEC**
- **Low Voltage Directive 73/23/EEC**

We hereby certify that the product specified below is in accordance with the provision of EU Machinery Directive 89/392/EEC - Annex II A, EU Electromagnetic Compatibility Directive 89/336/EEC and EU Low Voltage Directive 73/23/EEC.

The guidelines, harmonised standards, national standards and specifications which have been applied are listed below.

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Declaration of Contamination of Vacuum Equipment and Components

The repair and/or service of vacuum components will only be carried out if a correctly completed declaration has been submitted. Non-completion will result in delay.

The manufacturer could refuse to accept any equipment without a declaration.

This declaration can only be completed and signed by authorised and qualified staff:

1. Description of component:

- Equipment type/model: _____
- Code No.: _____
- Serial No.: _____
- Invoice No.: _____
- Delivery Date: _____

2. Reason for return:

3. Equipment condition

- Has the equipment been used?
yes no
- What type of pump oil was used?

- Is the equipment free from potentially harmful substances?
yes (go to section 5)
no (go to section 4)

4. Process related contamination of equipment

- toxic yes no
- corrosive yes no
- microbiological hazard*) yes no
- explosive*) yes no
- radioactive*) yes no
- other harmful substances yes no

*) We will not accept delivery of any equipment that has been radioactively or microbiologically contaminated without written evidence of decontamination!

Please list all substances, gases and by-products which may have come into contact with the equipment:

| Tradename Product name Manufacturer | Chemical name (or Symbol) | Danger class | Precautions associated with substance | Action if spillage or human contact |
|---|------------------------------|--------------|--|--|
| 1. | | | | |
| 2. | | | | |
| 3. | | | | |
| 4. | | | | |
| 5. | | | | |

5. Legally Binding Declaration

I hereby declare that the information supplied on this form is complete and accurate. The despatch of equipment will be in accordance with the appropriate regulations covering Packaging, Transportation and Labelling of Dangerous Substances.

Name of Organisation: _____

Address: _____ Post code: _____

Tel.: _____

Fax: _____ Telex: _____

Name: _____

Job title: _____

Date: _____ Company stamp: _____

Legally binding signature: _____

Zentrale/Headquarters

Pfeiffer Vacuum GmbH
Emmeliusstrasse 33
D-35614 Asslar
Telefon 06441/802-0
Telefax 06441/802-202
Hotline 06441/802-333
Internet:
<http://www.pfeiffer-vacuum.de>

Argentina

ARO S.A., Casilla de Correo 4890,
1000 Buenos Aires, telephone +54 / 1 331 3918,
telefax +54 / 1 331 3572

Australia

Balzers Australia Pty. Ltd., Level 1,
3, Northcliff Street, Milsons Point, NSW 2061,
telephone +61 / 2 9954 1925, telefax +61 / 2 9954 1939

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Great Britain

Pfeiffer Vacuum Ltd.
2-4 Cromwell Business Centre
Howard Way, Interchange Park
Newport Pagnell, MK16 9QS, United Kingdom
telephone +44 / 1 908 500615
telefax +44 / 1 908 500616

Greece

Analytical Instruments S.A., 1 Mantzarou St.,
GR-15451 Athens,
telephone +30 / 1 674 8973
telefax +30 / 1 674 8978

India

Pfeiffer Vacuum India Ltd.
25-E Nicholson Road, Tarbund
Secunderabad 500 009,
telephone +91 / 40 775 0014, telefax +91 / 40 775 7774

Israel

ODEM Scientific Applications,
2 Bergman Street, Science based park
Rehovod
telephone +972 8 9366 101,
telefax +972 8 9366 102

Italy

Pfeiffer Vacuum Italia S.p.a.
Via San Martino, 44 I-20017 RHO (Milano)
telephone +39 / 2 93 99 051, telefax +39 / 2 93 99 05 33

Japan

Hakuto Co. Ltd., C.P.O. Box 25,
Vacuum & Scientific Instruments Division
Tokyo Central 100-91,
telephone +81 / 3 32 258 910
telefax +81 / 3 32 259 009

Korea

Pfeiffer Vacuum Korea Ltd., 3F Haein Building 453,
Dokok-Dong, Kang Nam-Ku, Seoul, 135-270
telephone +82 / 2 3461 0671/5
telefax +82 / 2 3461 0676

Netherlands

Pfeiffer Vacuum Nederland BV
Veldzigt 30a, NL-3454 PW De Meern,
telephone +31 / 30 6666050, telefax +31 / 30 6662794

Peru

Ing. E. Brammertz s.r.l., José Pardo 182,
Apartado 173, PE-18 Miraflores,
telephone +51 / 1 445 8178
telefax +51 / 1 445-1931

Poland

Softtrade Sp.z.o.o, ul. Malwowa 35,
PL-60-175 Poznan, telephone +48 / 61 8677 168,
telefax +48 / 61 8677 111

Portugal

Unilaser Lda, Taguspark
Núcleo Central, sala nº 268, Estrada Cacém-
Porto Salvo, P-2780 Oeiras
telephone +351 / 1 421 7733
telefax +351 / 1 421 7744

Singapore

APP Systems Services Pte. Ltd, 2 Corporation Road
06-14 Corporation Place, Singapore 618494,
telephone +65 / 268 2024, telefax +65 / 268 6621

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Scope of represented countries

Bahrain, Egypt, Iraq, Iran, Jordan, Kuwait, Lebanon,
Lybia, Oman, Pakistan, Saudi-Arabia, Sudan, Syria,
Turkey, United Arab Emirates, Yemen

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Secotec S.A., Apartado 3452, Caracas 1010-A,
AVI - Applied Vacuum Industries GmbH