

Osmometers

► Semi-Micro Osmometer K-7400



Determination of total osmolality of biological fluids and other aqueous liquids

Typical applications for the Semi-Micro Osmometer K-7400 range from routine testing of physiological fluids such as blood serum and urine in clinics* up to quality control of food products like milk or isotonic soft drinks.

Principle

The Semi-Micro Osmometer K-7400 determines the osmolality by measuring the freezing point depression of the liquid sample. Osmolality describes the concentration of osmotic effective particles in solutions, independent from type, composition or electrical charge. The osmolality refers to the mass of the solution, meaning 1 kg of pure water.

Advantages

- Simple handling
- 2 to 4 point calibration
- Short measurement time (2–3 min)
- Reasonable price
- Fully automated measurement
- Automatic storage of measurement values
- Air cooling
- Cooling process is thermistor controlled
- Cooling chamber will not freeze



Medical* applications

- Poisoning
- Monitoring kidney functions
- Diabetes Insipidus
- Shock (trauma)
- Clearance control
- · Lactation control
- Dialysis control
- Control of infusion solutions in pharmacies
- Control of iso, hyper, and hypotonic solutions
- Blood control and control of other transfusion supplies
- Control of osmotic pressure in cell culture media for biotechnology and genetic engineering *In vitro diagnostics approval required for EU and other countries is in preparation.

Operation

150 µl of an aqueous solution or a biological fluid are super cooled below the melting point. Then a stirring needle (vibrator) automatically initiates crystallization. The crystallization warmth causes a temperature rise up to the melting point (below 0°C). The temperature difference to 0°C measured is proportional to the number of free particles in the solution. The measurement value in mOsmol is shown in the display.

Technical data

Measurement accuracy error ≤ 1% Measurement precision RSD < 1%

Calibration 2 calibration points necessary – up to 4 calibration points accepted

Sample volume 150 µl Measuring time 2 minutes Guaranteed range 0-2000 mOsmol Possible range 0-3500 mOsmol

Communication RS-232 interface connects to PC or printer A0591

Power supply 90–260 V, 47–63 Hz, max. 70 VA **Dimensions K-7400** 160 x 182 x 340 mm (WxHxD) **Printer** 76.8 x 77.4 x 39.3 mm (WxHxD)

Weight K-7400 5.0 kg Printer 0.15 kg

Ordering information

Order no.	Osmometer	Order no.	Accessories
A3707	Semi-Micro Osmometer K-7400	A0591	Printer
A3707-1	Semi-Micro Osmometer K-7400 with printer	A0840	Measuring Head for glass vials
A3709	Semi-Micro Osmometer K-7400	A0840-1	Measuring Head for plastic vials
	with Measuring Head for plastic vials	Y1240	Calibration Solution 300 mOsmol/kg
A3709-1	Semi-Micro Osmometer K-7400 with printer and with Measuring Head for plastic vials	Y1241	Calibration Solution 400 mOsmol/kg
		Y1243	Calibration solution 1 000 mOsmol/kg
		Y1248	Calibration solution 2000 mOsmol/kg
		A0913	Measuring Vials (glass), 1000 pcs
Technical data are subject to change without notice.		A0720	Measuring Vials (plastic), 1000 pcs

Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.

Wissenschaftliche Gerätebau Dr. Ing. Herbert Knauer GmbH Hegauer Weg 38 14163 Berlin, Germany



Your local distributor:

Printed on paper awarded the EU Ecolabel

Phone: +49-(0)30-80 97 27-0 Telefax: +49-(0)30-8 01 50 10 E-Mail: info@knauer.net Internet: www.knauer.net