



Centrifuge Bundle Guide



**SPECIAL
BUNDLE
PRICES**

Research Applications

High-Speed Micro Centrifuge Bundles



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30130866	FC5515 Micro	1	15 200	19 165 x g	1.5 / 2.0 mL	24
30642371	Angle Rotor 24 x 1.5/2.0 BIOSEALS	1				

A

List price: £2,006

BUNDLE DEAL: £1,595

Refrigerated



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30130868	FC5515R Micro Refrigerated	1	15 200	19 165 x g	1.5 / 2.0 mL	24
30642371	Angle Rotor 24 x 1.5/2.0 BIOSEALS	1				

B

List price: £4,919

BUNDLE DEAL: £3,995

For more information, or your own custom bundle, please contact the sales team at **Richmond Scientific** on 01257 270433

Haematocrit Bundle



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30393187	FC5513 Micro	1	13 500	16 085 x g	Capillary Tubes	24 Capillaries
30472304	Haematocrit Rotor 24 x Capillaries	1				

C

List price: £2,006

BUNDLE DEAL: £1,595

Microplate Bundle



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30314810	FC5714 Multi	1	4 500	2 716 x g	Plates	6 Standard Microplates or 2 Deepwell Plates
30314824	Swing Out Rotor 2x3 MTP with bucket	1				

D

List price: £4,145

BUNDLE DEAL: £3,295

For more information, or your own custom bundle, please contact the sales team at **Richmond Scientific** on 01257 270433

Cell Culture Bundles

Low Volume Cell Culture



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30314812	FC5718 Multi	1	5 000	3 885 x g		
30314823	Swing Out Rotor 4 x 100 mL Sealable	1			100 mL	4
30314882	Buckets – 2 x 50 mL Conical (2 pack)	2			50 mL	8
30314881	Buckets – 3 x 15 mL Conical (2 pack)	2			15 mL	12

E

List price: £5,669

BUNDLE DEAL: £4,549



Refrigerated

P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30314814	FC5718R Multi Refrigerated	1	5 000	3 885 x g		
30314823	Swing Out Rotor 4 x 100 mL Sealable	1			100 mL	4
30314882	Buckets – 2 x 50 mL Conical (2 pack)	2			50 mL	8
30314881	Buckets – 3 x 15 mL Conical (2 pack)	2			15 mL	12

F

List price: £8,779

BUNDLE DEAL: £6,995

For more information, or your own custom bundle, please contact the sales team at **Richmond Scientific** on 01257 270433

Cell Culture Bundles

Mid-Volume Cell Culture



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30314816	FC5816 Multi	1	4 500	37 355 x g		
30314828	Swing Out Rotor 4 x 250 mL	1			100 mL	4
30304375	Rack – 4 x 50 mL Conical (2 pack)	2			50 mL	16
30314583	Rack – 10 x 15 mL Conical (2 pack)	2			15 mL	40

G

List price: £7,009

BUNDLE DEAL: £5,495



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes per rotor
30314818	FC5816 Multi Refrigerated	1	4 500	37 355 x g		
30314828	Swing Out Rotor 4 x 250 mL	1			100 mL	4
30304375	Rack – 4 x 50 mL Conical (2 pack)	2			50 mL	16
30314583	Rack – 10 x 15 mL Conical (2 pack)	2			15 mL	40

H

List price: £10,359

BUNDLE DEAL: £8,295

Refrigerated

For more information, or your own custom bundle, please contact the sales team at **Richmond Scientific** on 01257 270433

Cell Culture Bundles

High-Volume Cell Culture



P/N	Description	Qty	RPM	RCF	Sample Size	Tubes /rotor
30553036	FC5916 Multi	1	4 000		750 mL	4
30553085	Swing Out Rotor 4 x 750 mL Sealable	1				
30602502	Bucket – 750 mL Flat Bottom (2 pack)	2				
30553130	Rack – 7 x 50 mL Conical (2 pack)	2				
30553134	Rack – 17 x 15 mL Conical (2 pack)	2				

I

List price: £9,410 **BUNDLE DEAL: £8,995**



Refrigerated

P/N	Description	Qty	RPM	RCF	Sample Size	Tubes /rotor
30553101	FC5916R Multi Refrigerated	1	4 000		750 mL	4
30553085	Swing Out Rotor 4 x 750 mL Sealable	1				
30602502	Bucket – 750 mL Flat Bottom (2 pack)	2				
30553130	Rack – 7 x 50 mL Conical (2 pack)	2				
30553134	Rack – 17 x 15 mL Conical (2 pack)	2				

J

List price: £12,190 **BUNDLE DEAL: £9,750**


For more information, or your own custom bundle, please contact the sales team at **Richmond Scientific** on 01257 270433



OHAUS Corporation

Headquartered in Parsippany, New Jersey, OHAUS Corporation is a global leader in manufacturing an extensive line of weighing products, laboratory equipment and analytical instruments that meet and surpass the weighing and measurement needs of a broad range of industries, including laboratory, industrial, educational, food preparation and retail markets. An ISO 9001:2015 manufacturer, OHAUS develops products that are precise, reliable, affordable and backed by industry-leading customer support.

For more information about any of the centrifuge bundles in this flyer, or to discuss custom pricing for your specific requirements, please contact the sales team at **Richmond Scientific** on 01257 270433



**Richmond
Scientific**

ENABLING
SCIENTIFIC DISCOVERY